
FOR FAECAL CONCENTRATION OF

HELMINTH OVA AND LARVAE / PROTOZOA CYSTS AND OOCYSTS



Health and Safety Benefits

- Totally enclosed/sealed process
- Reduced reagent volumes
- No cleaning required
- Single use, no sample contamination
- Ready to use systems available

Performance Benefits

- Optimum sample recovery
- Enhanced sample clarity
- Rapid four step process
- Human resources optimised
- Easy patient identification
- Fits all 50ml centrifuge buckets



PARASITOLOGY SINGLE USE IN VITRO DIAGNOSTIC DEVICE

APACOR

Midi Parasep[®] SI



Midi Parasep® S

STEP 1 - SAMPLE PREPARATION Fresh Samples

For empty Parasep[®], unscrew lid and add 8.0ml of fixative and one drop of surfactant (eg: Apacor Triton X solution) to the mixing chamber.

Alternatively use the reagent ready Midi Parasep[®] SF.

Introduce a level scoop of faecal sample to the fixative.


Mix in thoroughly with the Midi Parasep[®] SF spoon. If the sample is hard, break it up with the end of the spoon.


STEP 3 - CENTRIFUGATION

Invert the Midi Parasep[®] SF and centrifuge at 400g for two minutes (J. Clin. Microbiol. doi:10.1128/JCM.00838-15).

· · · ·

Midi Parasep® SF fits all 50ml centrifuge buckets.





STEP 2 - EMULSIFICATION

Procedure

Seal the Midi Parasep[®] SF by screwing in the filter/ sedimentation cone unit.

Vortex or shake to emulsify with the sedimentation cone pointing upwards.



STEP 4 - EXAMINATION Direct Method

Unscrew and discard the filter and mixing tube. Pour off all the liquid above the sediment.

Pipette one drop of sediment onto a slide, and cover with cover-slip. Alternatively, follow laboratory SOP for slide preparation.

OR

Semi-automated System - ParaSys™ Unscrew and discard the filter and mixing tube. Pour off all the liquid above the sediment.

Insert aspirator into the sediment, press DILUTE to add saline to the sediment. Shake or vortex to re-suspend sample. Press SAMPLE to draw 100µl into the ParaSlide™.



See label for storage conditions and expiry date. Please adhere to the following guidelines when handling Midi Parasep[®] SF. To avoid cross contamination the Midi Parasep[®] SF device should remain closed at all times except when introducing the sample or when retrieving the final concentrated sample for examination.

Midi Parasep[®] SF is available empty or reagent ready Please ask for details

Products can be ordered direct from Apacor or from an appointed distributor Visit our website for all the latest information www.apacor.com or email on: sales@apacor.com



UNIT 5 SAPPHIRE CENTRE FISHPONDS ROAD, WOKINGHAM BERKSHIRE, RG41 2QL, ENGLAND TEL: +44 (0)118 979 5566 FAX: +44 (0)118 979 5186









APA161 - V6.0 09/2017