**U.S. English**

*Product Number:*  
902093K

***Reveal<sup>®</sup> 3-D***  
*for Soy*

# Reveal<sup>®</sup> 3-D for Soy

Product Number: 902093K

## Introduction and Intended Use

Reveal<sup>®</sup> 3-D for Soy is uniquely designed with three lines of detection and can be used virtually anywhere to screen environmental swabs, rinses, food products and ingredients for the presence of significant levels of soy. It is not designed to be used quantitatively. The test's 3-D technology ensures greater reliability with screening than ever before.

Reveal 3-D for Soy is intended for use only in an industrial food manufacturing/preparation context. Because of the problems of adequately sampling and extracting soy, it is not suitable for the testing of foods to be consumed in the home or in a restaurant.

## Limit of Detection

The test is highly sensitive and specifically designed to screen for low parts per million, or ppm, levels of soy content in clean-in-place (CIP) rinses and environmental swabs and foods.

Using the environmental swabs supplied, levels of 2 µg/100 cm<sup>2</sup> of soy protein extract can be detected.

When analyzing rinses and foods, soy residues are detectable as low as 2.6 ppm soy protein (5 ppm total soy). The presence of cleaners and sanitizers in rinses, as well the type of food matrix being tested, can affect limit of detection.

## Test Performance

Reveal 3-D for Soy has undergone rigorous validation to evaluate the specificity, sensitivity, robustness and intra- and inter-batch variability of the test method on rinses, environmental swabs, and foods.

## Sample Compatibility

Reveal 3-D for Soy is designed to detect soy on environmental surfaces, in rinses and in foods. Although every effort has been made to validate as many matrices as possible, there are some that are not suitable for testing. **Note:** Reveal 3-D for Soy is not suitable for the detection of soy-derived lecithin.

Users should perform in-house matrix-specific spike and recovery validation work in conjunction with a validated laboratory assay, such as Veratox<sup>®</sup> for Soy to help confirm Reveal 3-D for Soy results. This process will highlight any problematic matrices encountered.

Existing commodity validations should be revalidated when suppliers or the manufacturing process has been changed. Please contact a NEOGEN representative or distributor for additional details.

## Materials Provided

Each Reveal 3-D for Soy kit contains the following:

1. An instruction leaflet
2. 1 sealed foil pouch, containing 10 beige Reveal 3-D Soy Test devices
3. 10 sachets containing Type 8 Extraction Buffer
4. 10 sample tubes and caps
5. 10 individually packaged, sterile swabs with break-off tips
6. 1 bottle of swab wetting solution

## Materials Required for Food Extraction But Not Provided

1. 3-D Food Buffer (NEOGEN item 8504)
2. 50 mL graduated cylinder tubes for extraction (NEOGEN item 9381)

## Sampling Technique

1. CIP rinse water  
As only a small sample of material is required for Reveal 3-D for Soy, it is important to test a representative portion of the liquid.
2. Environmental swabbing  
The swabs supplied are intended to be used for the collection of environmental samples from which the presence of soy can be tested. This method can be used to validate the adequacy of cleaning and/or to identify problem areas (e.g., unwanted buildup of soy in processing equipment).
3. Food samples  
As only a small sample of material is required for Reveal 3-D for Soy, it is important to test a representative portion of the sample.

## Precautions

1. For environmental testing in an industrial food manufacturing, preparation, or labeling enforcement context only.
2. Do not use any part of the test beyond the expiry date.
3. Do not open the foil pouch until just before use.
4. Ensure the foil bag is tightly sealed after removal of a device.
5. Always store the kit between 2–8°C (35–46°F). Avoid freezing.
6. Bring the kit to room temperature (18–30°C, 64–86°F) prior to use.

## Test Storage/Sample Stability

Store Reveal 3-D for Soy kit between 2–8°C (35–46°F) and use within the expiry date stated on the outer label. Extracted samples should be used within 3 hours of extraction.


## Limitations

A negative test from a surface or rinse cannot exclude the possibility that the food contains soy since it may be distributed unevenly on the surface or rinse and may be below the detection limit of the test with that specific sample.


Reveal 3-D for Soy is a qualitative test and should only be used as a preliminary screen for the presence of soy. The validity of results obtained with the test should preferably be viewed in conjunction with data from a quantitative assay, such as Veratox for Soy allergen.

Hydrolyzed and fermented proteins may not be detectable with antibody-based test methods for allergen testing. Although the proteins may be undetectable in the assay because of their properties, there still could be allergenic protein residue present.

Validate detection by testing a positive control from the source of the allergen or ingredient that is available in your manufacturing environment and that may pose an allergen contamination risk to ensure the test can detect the allergen of concern. Contact NEOGEN Technical Services for additional information.

1


3


## Rinse Sampling

For each sample to be tested, remove the following from refrigeration and allow to equilibrate at room temperature before use (20–30 minutes):


1. 1 sample tube
  2. 1 Type 8 Extraction Buffer
  3. 1 Reveal 3-D for Soy device (in foil pouch)
1. Carefully tear/cut off the top of the Type 8 Extraction Buffer and add the entire contents to the sample tube.
  2. Add 0.25 mL (250 ul) of sample to the sample tube.
  3. Secure the white cap and shake for one minute.


1


2

3

4

5

## Swab Sampling


For each sample to be tested, remove the following from refrigeration and allow to equilibrate at room temperature before use (20–30 minutes):

- 1 bottle of swab wetting solution
  - 1 sample tube
  - 1 Type 8 Extraction Buffer
  - 1 Reveal 3-D for Soy device (in foil pouch)
  - 1 sterile swab
1. Carefully tear/cut the Type 8 Extraction Buffer and add the entire contents to the sample tube.
  2. Estimate a swabbing area of approximately 10 cm x 10 cm. Alternatively, use the swab to collect samples of contamination from problem areas (e.g., processing equipment, filler heads).
  3. Gather the sample with the swab, using one of the following methods:


For dry surfaces: Remove a sterile swab from the packaging and wet with 2 drops of swab wetting material. Swab a 10 x 10 cm area by using a crosshatch technique, revolving the swab on the surface. Repeat this swabbing procedure using movements at right angles to those used in the first swabbing.

For wet surfaces: Remove a sterile swab from the packaging and swab a 10 x 10 cm area by using a crosshatch technique revolving the swab on the surface. Repeat this swabbing procedure using movements at right angles to those used in the first swabbing. Do not moisten swab prior to use.

4. Place the swab to the Extraction Buffer in the sample tube and carefully break off the moistened end at the pre-scored mark so that it remains in the tube.
5. Secure the cap of the sample tube, taking care to ensure that the stem does not prevent the tube from being properly sealed. Shake for 1 minute.

1


2


## Swab and Rinse Testing

1. Remove the lid and fill it with the liquid from the tube. Any froth should remain in the tube.
2. Dip the head of the Reveal 3-D device into the liquid in the lid. Ensure that the cavity is saturated with the liquid.
3. Leave the cavity saturated until liquid is seen running in the test window.
4. Place device on a flat surface and allow test to develop for 5 minutes.


**Note:** It is essential to place the device flat on a level surface as soon as the liquid has entered the test window to stimulate flow through the device. Additionally, the devices are pre-stripped with pale green loading dye in positions T (test), O (overload) and C (control). The loading dye assists with quality and manufacturing checks and does not impact test performance. The loading dye is removed from the test window as the sample flows through the device.

1

2

3

## Food Extraction and Testing

The extraction of food samples requires the use of NEOGEN's 3-D Food Buffer (NEOGEN item 8504).





From the test kit, remove the appropriate number of the following and allow to equilibrate at room temperature before use (20–30 minutes out of refrigerator):

3-D Food Buffer

Reveal 3-D for Soy device (in foil pouch)

1. Add 20 mL of 3-D Food Buffer into a new 50 mL tube (18 mL for liquid samples) or a similar container.
2. Weigh out 2 g (mL) of sample and add to tube with buffer.
3. Secure the cap and invert to suspend sample. Shake vigorously or vortex for 1 minute.
4. Tilt the tube at approximately a 45 degree angle. Dip the head of the Reveal 3-D device into the liquid. Ensure the cavity is saturated with the liquid, and that no liquid gets on the test window.
5. Leave the cavity saturated until liquid is seen running in the test window.
6. Place device on a flat surface and allow test to develop for 5 minutes.

**Note:** It is essential to place the device flat on a level surface as soon as the liquid has entered the test window to stimulate flow through the device. Additionally, the devices are pre-stripped with pale green loading dye in positions T (test), O (overload), and C (control). The loading dye assists with quality and manufacturing checks and does not impact test performance. The loading dye is removed from the test window as the sample flows through the device.

## Reading Results

Read the result 5 minutes after dipping. Observations after

6 minutes may be inaccurate due to overdevelopment of the device.

1. Negative result  
No line at position T (test): Level of soy is below the detection limit. (See Limitations section).
2. Positive result  
Any intensity of a red line at position T (test): Level of soy above detection limit.
3. High positive results  
No line is visible at position O (overload) and a red line is faintly visible or absent at position T: Sample is overloaded with soy.
4. Invalid results  
If no red line appears at position C (control), then the test may be invalid.

## Customer Service

NEOGEN Customer and Technical Services can be contacted through [NEOGEN.com](https://www.neogen.com) and product training is available by request.

## SDS Information Available

Safety data sheets are available for all test kits at [NEOGEN.com](https://www.neogen.com) or by calling 800.234.5333 or 517.372.9200.

## Terms and Conditions

NEOGEN's full terms and conditions are available [online](#).

## Warranty

NEOGEN makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, NEOGEN will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. NEOGEN shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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