

Aspergillus Galactomannan

Lateral Flow Assay

- 30 minute test time
- Simplified procedure
- High sensitivity and specificity
- Serum & BAL specimens
- Test more frequently /
Eliminate batch testing



Aspergillus GM Lateral Flow Assay

the Problem

Over 30 million people are at risk of invasive aspergillosis each year because of corticosteroid use or other therapies, and over 300,000 patients develop it annually.

The disease is common in high risk patients with:

- Hematological malignancies
- Allogeneic HSCT (stem cell transplant)
- Chemotherapy-induced neutropenia
- Solid organ transplant (primarily lung)

real Value

- Delivers accurate results in **30 minutes**, reducing time to proper treatment for patients.
- Test at risk patients more frequently.
- Simplified procedure with minimal specimen preparation.
- Eliminate delays in diagnostics due to batch testing

simple Procedure

SPECIMEN PREPARATION

Obtain 2 test tubes for each specimen: 1 screw cap, heat resistant centrifuge tube for the dilution
1 flat-bottom tube for running the test

1



Transfer 300 µl specimen in the screw cap, heat resistant centrifuge tube 1

2



Add 100 µl Sample Pretreatment Buffer to tube 1 (vortex as needed)

3




Place tube 1 on heat block for 6-8 minutes at 120°C

4



Centrifuge tube 1 at 10,000 - 14,000 x g for 5 minutes

5



Transfer 80 µL from tube 1 to tube 2

6



Add 40 µL of Aspergillus GM LFA Running Buffer to tube 2

7



Insert strip (↓ down) Wait for 30 min.

8



Read Test
1 line = negative
2 lines = positive

true Performance

The sōna AGM LFA displays high performance in both specimen types.

Serum	Asp Ag EIA	
	Pos	Neg
Asp GM LFA Assay	26	1
	6	116

Serum	Calculated	95% CI
% Agreement Pos	81%	64% - 93%
% Agreement Neg	99%	95% - 99.9%

BAL	Asp Ag EIA	
	Pos	Neg
Asp GM LFA Assay	25	3
	3	48

BAL	Calculated	95% CI
% Agreement Pos	89%	72% - 98%
% Agreement Neg	94%	84% - 99%