
 **NG·BIOTECH®**
ZA Courbouton, secteur 1
35480 GUIPRY - France
Tel : +33(0)2 23 30 17 83
Fax : +33(0)9 71 70 53 10
contact@ngbiotech.com
www.ngtest-covid-19.com

ENO228COV/Rév. 210107

NG·TEST®/SARS-CoV-2 Ag

Cassette

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4 / World Health Organization (WHO). Statement Regarding Cluster of Pneumonia Cases in Wuhan, China. Beijing: WHO; 9 Jan 2020. (Accessed 26 Jan 2020). <https://www.who.int/china/news/detail/09-01-2020-who-statement-regarding-cluster-of-pneumonia-cases-in-wuhan-china>.

5 / World Health Organization (WHO). Coronavirus. <https://www.who.int/health-topics/coronavirus>.

6 / Monita R, Darin C. Performance of Oropharyngeal swab testing compared with nasopharyngeal swab testing for Diagnosis of Coronavirus disease 2019 – United States, January 2020 – February 2020.

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INTENDED PURPOSE

The NG-Test® SARS-CoV-2 Ag is an immunochromatographic assay for the qualitative detection of SARS-CoV-2 antigen in nasopharyngeal and oropharyngeal specimens. The NG-Test® SARS-CoV-2 Ag is for professional use only and is intended for use as an aid in the diagnosis of SARS-CoV-2 infection. NG Biotech determined the performance characteristics of the assay based on nasopharyngeal specimens. Further, the assay is compatible with oropharyngeal specimens. Note that the concentration of viral antigen may vary in the course of the disease and might fall below the detection threshold of the test, particularly when the duration of disease symptoms increases (specimens collected within 7 days of the first symptom are generally more likely to give a positive result). Infectivity cannot be excluded on the basis of negative test results alone.

SUMMARY

Early in January 2020, a new coronavirus (named "SARS-CoV-2") was identified as being the infectious agent causing a viral pneumonia outbreak in Wuhan, China, where the first cases experienced symptoms in December 2019. The disease it causes was named "coronavirus disease 2019" (in short "COVID-19"). Coronaviruses are enveloped RNA viruses widely present among human beings, other mammals, and birds, and lead to respiratory, enteric, hepatic, and neurological diseases. Six coronavirus species are known to cause human diseases. Four viruses - 229E, OC43, NL63, and HKU1 - are common and typically cause common cold symptoms in immunocompetent people. The two other strains - severe acute respiratory syndrome coronavirus (SARS-CoV) and Middle East respiratory syndrome coronavirus (MERS-CoV) - are of zoonotic origin and have been associated with some potentially life-threatening diseases. Common COVID-19 infection signs include respiratory symptoms, fever, cough, shortness of breath, and breathing difficulties. In more serious cases, the infection may cause pneumonia, severe acute respiratory syndrome, renal failure, and even death. The standard recommendations to prevent infection spread include regular hand washing, cover of mouth and nose when coughing and sneezing. Avoid any close contact with any person experiencing respiratory disease symptoms such as cough and sneezing.

PRINCIPLE

During the test, the sample reacts with colloidal gold nanoparticles conjugated to monoclonal antibodies directed against SARS-CoV-2 nucleocapsid protein. Then, the mixture migrates through the membrane by capillary action, and reacts with another monoclonal antibody directed against SARS-CoV-2 nucleocapsid protein printed at the test line (T) level. If the SARS-CoV-2 concentration present in the sample is higher than the detection limit, a coloured line appears on line T.

A colour line should always be present on the control line marked "C", indicating sufficient amount was applied and the test worked properly.

REAGENTS AND MATERIALS SUPPLIED

Each kit contains:

- 20 cassettes individually packaged in aluminum pouches with desiccant
- 20 extraction tubes and 20 caps *
- 2 extraction buffers in plastic dropper bottles
- 20 sterile swabs*
- 1 tube holder
- 1 package insert

*Due to potential shortages of ancillary medical products related to COVID-19, the nasopharyngeal swabs and extraction tubes with caps provided in the kit may be subject to change.

MATERIALS REQUIRED BUT NOT SUPPLIED

Timer • Personal protective equipment • Micropipette • Microtube

PRECAUTIONS

- For professional in vitro diagnostic use only.
- It is recommended to use the test under a type-II microbiological safety cabinet to protect the person handling SARS-CoV-2 virus.
- Do not use the test if the pouch is torn or damaged.
- If the pouch was stored between 4 and 8°C, wait at least 10 minutes so that the test reaches room temperature.
- The test cassette should remain in the sealed pouch until use.
- Once the aluminium pouch is open, carry out the test rapidly.
- Do not reuse the device, the swab, or the extraction tube.
- Handle all the samples as if they contained infectious agents. Follow the established precautions against microbiological risks throughout all the procedures, and follow the standard procedures to collect and dispose of the samples appropriately.
- Wear protective clothing such as masks, lab coats, disposable gloves, and protective eyewear when analysing the samples.
- Make sure an appropriate amount of sample is used to carry out the tests. If the sample volume is too high or too low, the results may deviate.
- The used test should be disposed of in accordance with local regulations.
- Do not eat, drink, or smoke within the area where samples and kits are handled.

STORAGE AND STABILITY

Store the test in the sealed aluminium pouch between 4 and 30°C. The test is stable until the expiry date printed on the pouch. The test device should remain in the sealed pouch until use. Do not freeze. Do not use after the expiry date.

SPECIMEN COLLECTION

- 1 / Wear protective gloves, eyewear, and mask.
- 2 / Bring tests to room temperature prior testing.
- 3 / Place the clean extraction tube labelled with the patient ID into the provided tube holder.
- 4 / Holding the buffer bottle vertically and without touching the edge of the tube, add 10 drops into the extraction tube (~300µL). In case of foam or bubbles, use the second buffer bottle.
- 5 / Collect the nasopharyngeal or oropharyngeal specimen using the standard laboratory method and the swab provided in the kit.

Use the swab provided in the kit, **do not use any other swabs**.

**NASOPHARYNGEAL SPECIMEN**

- 1 / Keep the patient's head slightly tilted backwards.
- 2 / Insert the swab into the nostril and gently push it as far as possible along the floor of the nasal cavity, parallel to the palate, as described in the picture:
 - Swab the epithelium of the nasal cavity with rotating movements for 15 seconds, scraping the walls of the nasopharynx to collect cells.
 - Slowly remove the swab and place it in the extraction tube.

Notes :

Swabbing can be performed in both nostrils of the same patient using the same swab. However, this is not necessary if the swab is properly saturated with fluid from the first sampled nostril.
Do not use samples that are contaminated with blood.

**OROPHARYNGEAL SPECIMEN**

- 1 / Insert the swab into the posterior pharynx and gently brush it several times against the reddened posterior pharyngeal wall and both tonsillar pillars.
- 2 / Slowly remove the swab and place it in the extraction tube.


Notes :

Literature suggests that nasopharyngeal swabs may more accurately detect the amount of SARS-CoV-2. Although oropharyngeal samples can be used with this kit, it is recommended to use nasopharyngeal samples if testing environment allows.

**POINT OF CARE PROCEDURE**

The following procedure is intended to collect fresh samples directly in patients according to the above mentioned procedure.

- 1 / Unload the swab immediately after sampling by rotating it in the buffer 5-10 times and squeezing it to extract as much volume of sample as possible.
- 2 / Before removing the swab, squeeze it against the tube walls to release as much liquid as possible.
- 3 / Put the swab back into its original package and dispose of it according to local regulations.
- 4 / Remove the cassette from the aluminium pouch. Place the test on a clean and flat surface.
- 5 / Immediately after extraction, attach a dropper cap to the extraction tube, add 3 drops (120µL) of the mixture in the "S/R" well of the cassette, and then start the timer.
- 6 / Read the result after exactly 15 minutes. Do not read after more than 20 minutes to avoid false positive results.
- 7 / Dispose of the test and the tube in accordance with local regulations.

WAIT TO GET THE RESULT 15 min.

PROCEDURE IN A LABORATORY

The procedure below is intended to use samples that are already collected and diluted in viral transport medium. The viral transport medium (VTM), without denaturing agents, can be used to store samples for subsequent detection of the antigen.

The following media are recommended: UTM (Copan, ref 346C), saline solution, deionised water.

Use the recommended VTM_s, **and do not use other VTM_s**.

- 1 / Wear protective gloves, eyewear, and mask.
- 2 / Bring the test to room temperature prior testing.
- 3 / Add 4 drops buffer in a microtube**.
- 4 / With a micropipette**, transfer 120µL VTM-stored nasopharyngeal sample in the extraction tube.
- 5 / Homogenise the mixture through repeated aspirations with the micropipette.
- 6 / Transfer 120µL mixture in the cassette "S/R" well, and start the timer.
- 7 / Read the result after exactly **15 minutes**. Do not read after more than 20 minutes to avoid false positive results.
- 8 / Dispose of the test and the tube according to local regulations.

**Accessories not provided in the kit.


RESULT INTERPRETATION

NEGATIVE

Negative result

Only one red line appears within the control area (C). The sample does not contain SARS-CoV-2 or the virus concentration cannot be detected by the test. This should be interpreted as a negative result.

POSITIVE

Positive result

Two red lines appear, one in the control area (C) and one in the test area (T). The sample contains SARS-CoV-2 virus and should be interpreted as a positive result.

NOTE: The intensity of the red test line (T) may vary depending on SARS-CoV-2 virus level in the sample. A low-intensity line should be considered as a positive result.

INVALID

Invalid result

If the control line (C) does not appear, the test result is invalid. Most often, the control line does not appear because of insufficient sample volume or an incorrect procedure. The kit may have been deteriorated. Repeat the procedure using a new test. If the problem persists, do not reuse the kit and contact your distributor.

QUALITY CONTROL

An internal control is included in the test. When the control line appears, it confirms the sample volume was sufficient and the procedure was correctly applied. No control sample is provided with this kit.

LIMITATIONS

- 1 / NG-Test® SARS-CoV-2 Ag s for professional *in vitro* diagnostic use only.
- 2 / This test does not allow to establish the quantitative SARS-CoV-2 value in the sample.
- 3 / This test allows to confirm the presence of SARS-CoV-2 nucleocapsid antigen in a symptomatic patient. In case of negative result **it shouldn't be used as the only diagnostic criterion for SARS-CoV-2 infection**.
- 4 / As for any diagnostic test, the doctor should interpret the results taking into account the other clinical data.
- 5 / If the test result is negative, a confirmation test by PCR should be performed.
- 6 / The test will produce a negative result if the concentration of SARS-CoV-2 virus in the sample is lower than the limit of detection or if the nasopharyngeal specimen collection was not performed adequately.
- 7 / The detection limit of the immunodassay is the same among all qualified specimen e.g. nasopharyngeal and oropharyngeal samples. Nevertheless, professional preanalytical sample collection is essential for optimum results.

PERFORMANCES

1 / Evaluation of NG-Test® SARS-CoV-2 Ag

NG-Test® SARS-CoV-2 Ag was assessed with 200 nasopharyngeal samples confirmed by RT-PCR, including 100 positive samples and 100 negative samples. The overall specificity was evaluated at 100% with IC95% [96-100%].

The results are presented in th following tables:

Ct* VALUE	RT-PCR			-
	POSITIVE	NEGATIVE	OVERALL	
POSITIVE	75	81	83	0
NEGATIVE	4	12	17	100
OVERALL	79	93	100	100
SENSITIVITY			SPECIFICITY	
95%	87%	83%	100%	
IC 95% [88–99%]	IC 95% [79–93%]	IC 95% [74–90%]	IC 95% [96–100%]	

*Ct = Cycle threshold

2 / Limit of detection

The limit of detection was evaluated at 1ng/mL with recombinant nucleocapsid protein. Higher concentration will lead to positive results when processed on the assay.

3 / Repeatability and reproducibility

Repeatability and reproducibility of NG-Test® SARS-CoV-2 Ag were evaluated with positive an negative samples. No significant difference was observed with between lots, between operators and inter series results.

4 / Interferences

The following transport media have no impact on NG-Test® SARS-CoV-2 Ag : UTM (Copan, ref 346C), physiological serum, demineralized water.

5 / Cross reactions

a. Human corona viruses

No cross-reaction was observed with the viral nucleoprotein antigens of the coronavirus HCoV-NG3, MERS-CoV, HCoV-OC43, HCoV-229E strains up to a concentration of 1µg/mL (i.e. equivalent to about 12.04 x 10⁹ particles/mL). The viral nucleoprotein of SARS-CoV coronavirus strain cross-reacts with the NG-Test® SARS-CoV-2 for concentrations superior or equal to 0.5 ng/mL (i.e. superior or equal to about 6.02 x 10⁶ particles/mL).

b. UK new strain (VOC202012/01)

The results from the study results conducted on the recombinant nucleocapsid protein from VOC202012/1 variant (also known as B.1.1.7 or „UK strain“) show that the monoclonal antibodies used in the NG Test® SARS-CoV-2 Ag are able to detect this new variant of concern.

c. Miscellaneous

Cross reactions on NG-Test® SARS-CoV-2 Ag were evaluated by testing the following bacteria and viruses: enterovirus, respiratory syncytial virus, adenovirus, influenza B, parainfluenza, mycoplasma pneumoniae, Bordetella pertussis, human metapneumovirus, influenza A H1 et influenza A H3, influenza A H1 2009. Those viruses have not impacted on the results of the NG-Test® SARS-CoV-2 Ag.

6 / Matrix effect

A matrix equivalency study has been performed with negative nasopharyngeal and oropharyngeal specimens spiked with inactivated virus and recombinant protein and shows equivalent results.

BIBLIOGRAPHY

- 1 / Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses. *Nat. Rev Microbiol* 2019;17:181-192. PMID:30531947 DOI:10.1038/s41579-018-0118-9
- 2 / Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses. *Trends Microbiol* 2016;24:490-502. PMID:27012512 DOI:10.1016/j.tim.2016.03.003.
- 3 / Weiss SR, Leibowitz JL. Coronavirus pathogenesis. *Adv Virus Res* 2011;81:85-164. PMID:22094080 DOI:10.1016/B978-0-12-385885-6.00009-2.
- 4 / World Health Organization (WHO). Statement Regarding Cluster of Pneumonia Cases in Wuhan, China. Beijing: WHO; 9 Jan 2020. (Accessed 26 Jan 2020). <https://www.who.int/china/news/detail/09-01-2020-who-statement-regarding-cluster-of-pneumonia-cases-in-wuhan-china>.
- 5 / World Health Organization (WHO). Coronavirus. <https://www.who.int/health-topics/coronavirus>.
- 6 / Monita R, Darin C. Performance of Oropharyngeal swab testing compared with nasopharyngeal swab testing for Diagnosis of Coronavirus disease 2019 – United States, January 2020 – February 2020.

SYMBOLS

	Sufficient content for 20 tests		Read the instructions before use		Expiry date		Temperature limits (4°C-30°C)
	In vitro Diagnostic Medical Device		Manufacturer		Do not reuse		Do not use if package is damaged
	Batch number		Keep dry		Product reference		

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