● CHROMagar™
C.difficile


For detection of *Clostridium difficile*

Plate Reading

- *C. difficile*
→ colourless and fluorescent under UV light at 365 nm
- Other bacteria
→ colourless, no fluorescent or inhibited

For detection of *Clostridium difficile*

Background

Clostridium difficile (*C. difficile*) is the leading cause of nosocomial infectious diarrhea in adults. These infections occur mostly in patients who have both medical care and antibiotic treatment. The symptoms of *C. difficile* infection are fever, abdominal cramps and severe diarrhea leading to death. In the United States, nearly 250,000 people each year develop *C. difficile* infections with at least 14 000 deaths (CDC estimate, 2013). Due to the emergence of highly toxigenic *C. difficile* strains, these infections have become more frequent and more difficult to treat in the last years.

Although PCR has become the leading *C. difficile* detection technique, culture is essential for strain typing and antimicrobial susceptibility testing. CHROMagar™ *C. difficile* is a new fluorogenic culture medium, extremely sensitive and selective, especially designed to simplify and speed up (24h) the culture of *C. difficile*.

Medium Performance

- 1 RAPID DETECTION COMPARED TO TRADITIONAL MEDIA**
big colonies (around 2mm) of *C. difficile* after only 24h of incubation in anaerobic atmosphere, contrary to traditional media requiring 48h.
- 2 HIGH SENSITIVITY AND SPECIFICITY**
C. difficile is detected by characteristic fluorescent colonies (under UV light at 365 nm) and the specimen's flora largely inhibited
- 3 POLYVALENCE**
This medium can be used for clinical specimens as well as environmental samples.

Specificity / Sensitivity: ≈ 100%*

* Sensitivity from scientific studies: Gaillot O. et al. 100% (40/40), Van Broek et al. 100% (95/95)

Medium Description

Powder Base	Total	54.7 g/L
	Agar	15.0
Supplement (Included in the pack)	Peptone & Yeast extract	25.0
	Salts	9.0
	Growth factors	4.0
	Chromogenic mix.....	1.7
	Storage at 15/30°C - pH: 7.8 ± 0.2	
	Shelf Life	2 years
	Powder form.....	325mg/L
Storage at 2/8°C	Shelf Life	2 years

Usual Samples	Stool, Environmental
Procedure	Direct streaking or after an appropriate enrichment step of the sample. Incubation at 37°C for 24h. Anaerobic condition.

Scientific Publications on this product: available on www.CHROMagar.com
Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com

Quality Control Strains

<i>C. difficile</i> ATCC® 43255	colourless and fluorescent
<i>C. perfringens</i> ATCC® 13124	inhibited
<i>E. faecalis</i> ATCC® 29212	inhibited
<i>E. coli</i> ATCC® 25922	inhibited
<i>C. albicans</i> ATCC® 10231	inhibited

ATCC® is a registered trademark of the American Type Culture Collection

Order References

Please use these product references when contacting your local distributor:
5000 ml pack CD122
(included in this reference: powder base CD122(B) + supplement CD122(S))

Manufacturer: CHROMagar
4 place du 18 juin 1940 75006 Paris - France
Email: CHROMagar@CHROMagar.com
Website: www.CHROMagar.com
Find your nearest distributor on www.CHROMagar.com/contact

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