● CHROMagar™
ESBL

**For overnight detection
of Gram-negative bacteria producing
Extended Spectrum Beta-Lactamase**

For overnight detection of Gram-negative bacteria producing Extended Spectrum Beta-Lactamase Background

ESBL (Extended Spectrum β-Lactamases) are enzymes that mediate resistance to penicillins, extended-spectrum third generation cephalosporins (C3G) and monobactams. ESBL-producing *Enterobacteriaceae* started to appear in the 1980s, and have since emerged as some of the most significant hospital-acquired infections with *Escherichia coli* and *Klebsiella spp.* being the main actors, but other Gram-negative species have also been observed. Emergence of ESBL-producing isolates has important clinical and therapeutic implications:

- resistance determinants for ESBL production are carried on plasmids that can be easily spread from organism to organism.
- the spread of resistance toward extended-spectrum cephalosporins may lead to increased prescription of more broad-spectrum and expensive drugs.
- these resistant isolates may escape detection with routine susceptibility testing performed by a clinical microbiology laboratory, which can result in adverse therapeutic outcomes.

Therefore, the early detection of ESBL-producing bacteria carriers is important to minimise their impact and the spread of infections and customise therapeutic patient treatment.

Medium Performance

CHROMagar™ ESBL allows the detection of ESBL-producing bacteria while inhibiting the growth of other bacteria, including most of those carrying AmpC type resistance. This is an important feature because intrinsic AmpC resistance has less epidemic relevance, but often leads to ESBL false positive reading in the classical testing methods.

Detection of resistant isolates is difficult based on routine susceptibility testing performed by a clinical microbiology laboratory, while with CHROMagar ESBL:


Plate Reading

- *E.coli* ESBL → dark pink to reddish
- *Klebsiella, Enterobacter, Citrobacter* ESBL → metallic blue (+/- red halo)
- *Proteus* ESBL → brown halo
- *Pseudomonas* ESBL → translucent cream to blue
- *Acinetobacter* ESBL → cream, opaque
- Other → inhibited

1 FAST RESULTS

Detection after overnight incubation

2 SPECIES DIFFERENTIATION

thanks to the chromogenic performances of supplemented CHROMagar™ Orientation. Indeed, the product is composed of a powder base CHROMagar™ Orientation and a supplement to enhance ESBL-producing bacteria.

3 HIGH SENSITIVITY (99,2%*) HIGH SPECIFICITY (89%*)

*«Evaluation of a chromogenic medium for extended-spectrum beta-lactamase (ESBL)-producing *Enterobacteriaceae*» Philippe Lagacé-Wiens et al. University of Manitoba, Canada. ECCMID Poster 2010.

4 TIME AND WORKLOAD SAVINGS

Direct culture from specimen. There is no need of a selective pre-enrichment.

5 FLEXIBILITY

CHROMagar™ ESBL is supplied with a shelf-life of about 2 years. This allows flexibility of use, whether in an epidemic situation with many patients to screen, or whether for random surveillance of cultures.

Medium Description

Powder Base (CHROMagar Orientation)	Total	33 g/L
	Agar	15.0
	Peptone and yeast extract	17.0
	Chromogenic mix	1.0
	Storage at 15/30°C - pH: 7.0 +/-0.2	
	Shelf Life	2 years
+		
CHROMagar ESBL Supplement (included in the pack)	Selective mix (Powder form)	0.57 g/L
	Storage at 2/8°C	
	Shelf Life	2 years

Usual Samples	stools, urine
Procedure	Direct Streaking. Incubation 18-24h at 37°C. Aerobic conditions

Scientific Publications on this product: available on www.CHROMagar.com
Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com

Quality Control Strains

<i>E. coli</i> ESBL CIP 103982	reddish, small colonies
<i>K. pneumoniae</i> ATCC® 700603	metallic blue
<i>E. faecalis</i> ATCC® 29212	inhibited
<i>P. aeruginosa</i> ATCC® 10145	inhibited
<i>E. coli</i> ATCC® 25922	inhibited
<i>C. albicans</i> ATCC® 60193	inhibited
<i>S. aureus</i> ATCC® 25923	inhibited

ATCC® is a registered trademark of the American Type Culture Collection

Order References

Please use these product references when contacting your local distributor:

- 5000 ml pack ESRT2 (includes a base powder RT412 and a supplement ES372)
- 25 L pack ESRT3-25 (includes a base powder RT413-25 and a supplement ES373-25)
- Bulk on request

Manufacturer: CHROMagar
4 place du 18 juin 1940 75006 Paris - France
Email: CHROMagar@CHROMagar.com
Website: www.CHROMagar.com
Find your nearest distributor on
www.CHROMagar.com/contact

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