● CHROMagar™
KPC

**For detection of gram-negative bacteria
with a reduced susceptibility
to most of the carbapenem agents**



Plate Reading

-  • *E. coli* Carbapenem^R
→ dark pink to reddish
-  • *Klebsiella*, *Enterobacter*,
Citrobacter Carbapenem^R
→ Metallic blue (+/- red halo)
-  • *Pseudomonas* Carbapenem^R
→ translucent
cream to blue
- *Acinetobacter* Carbapenem^R
→ cream, opaque
- Carbapenem^S strains
→ inhibited

For detection of gram-negative bacteria with a reduced susceptibility to most of the carbapenem agents

Background

Worldwide reports of resistance to carbapenems found in enterobacteriaceae are a major health concern, specially in the case where the reducing susceptibility mechanism is the production of enzymes like KPC, OXA or MBL (ex: the recently reported NDM-1).

Carbapenems are the last resort in treating many serious gram-negative infections. However, production of these enzymes results in resistance to penicillins, cephalosporins (i.e., cefepime, ceftriaxone), carbapenems (i.e., meropenem, ertapenem), and aztreonam, thereby making these pathogens truly multidrug-resistant and making their treatment very challenging.

« KPC-producing bacteria have demonstrated a remarkable ability to disseminate with inter-facility, interstate, and international transmission having been documented. » CDC 2008-R-24. Thus, in order to limit the spread of these serious pathogens, rapid detection, followed by implementation of adequate infection control methods, is essential.

Medium Performance

- 1 DETECTION AFTER OVERNIGHT INCUBATION**
Detection of gram-negative bacteria expressing a reduced susceptibility to antibiotics of the carbapenem family.
- 2 TIME AND WORKLOAD SAVINGS**
There is no need for a selective pre-enrichment medium. Direct plating of the sample is possible.
- 3 SHORT INCUBATION**
Requires only 18-24 hours of incubation.
- 4 UNIQUE**
Unique commercially available culture medium for the direct isolation of bacteria harbouring reduced susceptibility to carbapenems
- 5 FLEXIBILITY**
CHROMagar™ KPC Supplement is supplied with a shelf-life of about 2 years. This allows for flexibility of use, whether in an epidemic situation with many patients to screen, or whether for random surveillance of cultures.

Medium Description

Powder Base CHROMagar Orientation	Total	33 g/L
	Agar	15.0
	Peptone and yeast extract	17.0
	Chromogenic mix	1.0
	Storage at 15/30°C - pH: 7.0 +/-0.2	
	Shelf Life	2 years
+		
CHROMagar KPC Supplement (included in the pack)	Selective mix (Powder form)	0.4 g/L
	Storage at 2/8°C	
	Shelf Life	2 years

Usual Samples	stools, urine
Procedure	Direct Streaking. Incubation 18-24h at 37°C Aerobic conditions.
Scientific Publications on this product: available on www.CHROMagar.com Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com	

Quality Control Strains

<i>K. pneumoniae</i> ATCC® BAA 1705.....	steel blue
<i>K. pneumoniae</i> NCTC 13438	steel blue
<i>E. faecalis</i> ATCC® 29212	inhibited
<i>K. pneumoniae</i> ATCC® 13883	inhibited
<i>S. aureus</i> ATCC® 25923	inhibited
<i>C. albicans</i> ATCC® 60193	inhibited

ATCC® is a registered trademark of the American Type Culture Collection

Order References

Please use these product references when contacting your local distributor:


- 5000 ml pack KPRT2 (includes a base powder RT412 and a supplement KP102)
- 25 L pack KPRT3-25 (includes a base powder RT413-25 and a supplement KP103-25)

Manufacturer: CHROMagar
4 place du 18 juin 1940 75006 Paris - France
Email: CHROMagar@CHROMagar.com
Website: www.CHROMagar.com
Find your nearest distributor on
www.CHROMagar.com/contact

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