CHROMagar™ COL-APSE

Instructions For Use
Available in several languages

NT-EXT-100

Version 4.0

ENGLISH

English Version

FRANCAIS


Version Française

ESPAÑOL

Version Español

DEUTSCH

Deutsch Version

MEDIUM PURPOSE

Chromogenic medium for detection of Colistin resistant Gram-negative bacteria.

Polymyxin E (colistin) and B are increasingly used as antimicrobials in the treatment of multi-drug resistant bacterial infections. Polymyxin resistance, although intrinsic in Gram-positive and some Gram-negative species (*Proteus*, *Morganella*, *Serratia*), is now a problem in a number of other pathogens (*Acinetobacter baumannii*, *Pseudomonas aeruginosa*, *Escherichia coli*, *Salmonella enterica*, *Klebsiella pneumoniae*).¹

Resistance arises due to mutations / insertions in genes involved in LPS biosynthesis (lpx, pmrA/B, mgrB, phoP/Q) and / or the acquisition of phosphoethanolamine transferases (PEtN). Of great concern is the recently described plasmid-encoded PEtN, MCR-1, now found worldwide in a range of animal, human and environmental bacterial isolates.¹

1- Novel Chromogenic Culture Media (CHROMagar™ COL-APSE) for the Isolation and Differentiation of Colistin Resistant Gram-negative Pathogens. ECCMID 2017.

COMPOSITION

The product is composed of a powder base (B) and a supplement (S).

Product	=	Base (B)	+	Supplement (S)
Total		42.5 g/L		2 ml/L
Composition		Agar 15.0 Peptones 20.0 Salt 5.0 Chromogenic and selective mix 0.8 Growth factors 1.7		Growth factors mix
Aspect		Powder Form		Liquid Form
STORAGE		15-30 °C		15-30 °C
FINAL MEDIA pH		7.1 +/- 0.2		

PREPARATION (Calculation for 1 L)

Step 1

Preparation of
Base + S

- Disperse slowly 42.5 g of powder base in 1 L of purified water.
- Add 2ml of CHROMagar™ COL-APSE supplement S into slurry.
- Stir until the agar is well thickened.
- Heat and bring to boiling (100 °C) while swirling or stirring regularly. DO NOT HEAT TO MORE THAN 100 °C. DO NOT AUTOCLAVE AT 121 °C.

Warning 1: If using an autoclave, do so without pressure.

Advice 1: For the 100 °C heating step, mixture may also be brought to a boil in a microwave oven: after initial boiling, remove from oven, stir gently, then return to oven for short repeated bursts of heating until complete fusion of the agar grains has taken place (large bubbles replacing foam).

- Cool in a water bath to 45-50 °C, swirling or stirring gently to homogenize.

Advice 2: in case of product samples containing a high load of *Proteus*, Cefixime can be added after the previous step at 0.05 mg/L.

Step 2

Pouring

- Pour into sterile Petri dishes.
- Let it solidify and dry.

Storage

- Store in the dark before use.
- Prepared media plates can be kept for one day at room temperature.
- Plates can be stored for up to 1 month under refrigeration (2/8 °C) if properly prepared and protected from light and dehydration.

INOCULATION

Related samples can be processed by direct streaking on the plate.

- If the agar plate has been refrigerated, allow to warm to room temperature before inoculation.
- Streak sample onto plate.
- Incubate in aerobic conditions at 37 °C for 18-24 hours.

Typical Samples


e.g. stools, rectal swabs,
veterinary (caecal, etc),
environmental

Direct streaking
or spreading technique

INTERPRETATION

Microorganism	Typical colony appearance
COL-R <i>E. coli</i>	→ dark pink to reddish
COL-R Coliforms	→ metallic blue
COL-R <i>Pseudomonas</i>	→ translucent, +/- natural pigmentation cream to green
COL-R <i>Acinetobacter</i>	→ cream, opaque
Other COL-R Gram negative	→ colourless, natural pigmentation
COL-S Gram negative bacteria	→ inhibited
Gram positive bacteria, yeasts	→ inhibited

Typical colony appearance

The pictures shown are not contractual.

PERFORMANCE & LIMITATIONS

- Species final identification may require additional testing such as biochemical tests.
- Resistance to colistin for *E. coli*, *Klebsiella*, *Enterobacter* and *Salmonella* can be confirmed by microdilution method.

QUALITY CONTROL

Please perform Quality Control according to the use of the medium and the local QC regulations and norms. Good preparation of the medium can be tested, isolating the reference strains below:

Microorganism	Typical colony appearance
MCR-1 <i>E. coli</i> NCTC 13846	→ dark pink to reddish
COL-S <i>E. coli</i> ATCC® 25922 (WDCM 00013)	→ inhibited
<i>E. faecalis</i> ATCC® 29212 (WDCM 00087)	→ inhibited

WARNINGS

- Do not use plates if they show any evidence of contamination or any sign of deterioration.
- Do not use the product beyond its expiry date or if product shows any evidence of contamination or any sign of deterioration.
- For Laboratory use. This laboratory product should be used only by trained personnel in compliance with good laboratory practices.
- Any change or modification in the procedure may affect the results.
- Any change or modification of the required storage temperature may affect the performance of the product.
- Inappropriate storage may affect the shelf life of the product.
- Recap the bottles/vials tightly after each preparation and keep them in a low humidity environment, protected from moisture and light.
- For a good microbial detection: collection and transport of specimen should be well handled and adapted to the particular specimen according to good laboratory practices.

DISPOSAL OF WASTE

After use, all plates and any other contaminated materials must be sterilized or disposed of by appropriate internal procedures and in accordance with local legislations. Plates can be destroyed by autoclaving at 121°C for at least 20 minutes.

REFERENCES

Please refer to our website page «Publications» for scientific publications about this particular product.
 Web link: <http://www.chromagar.com/publication.php>

IFU/LABEL INDEX

- Quantity of powder sufficient for X liters of media
- Expiry date
- Required storage temperature
- Store away from humidity

Need some Technical Documents?

- Available for download on www.CHROMagar.com
- Certificate of Analysis (CoA) --> One per Lot
- Material Safety Data Sheet (MSDS)

	Pack Size		Ordering References		Base (B)	Supplement (S)
	5000 ml		CO262	=	CO262(B) + Weight: 212.5 g	CO262(S) Volume: 10 ml

CHROMagar™ and Rambach™ are trademarks created by Dr A. Rambach
 ATCC® is a registered trademark of the American Type Culture Collection
 NT-EXT-100 V4.0 / EN 16-Jul-18

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