

Specification

Non-nutritive medium used for transporting of clinical samples and for virus preservation before analysis according to the CDC formulation.

Presentation

20 Tubes
Tubo PP
with: 3 ± 0,1 ml

Packaging Details

Polypropylene tubes, Non injectable cap. - 20 tubes per box.

Shelf Life

12 months

Storage

2-8 °C

Composition

Composition (g/l):

Sodium chloride.....	8.00
Potassium chloride.....	0.40
Monopotassium phosphate.....	0.06
Dextrose.....	1.00
Disodium phosphate.....	0.05
Magnesium sulfate.....	0.10
Calcium chloride.....	0.14
Sodium bicarbonate.....	0.35
Fetal bovine serum.....	20.0 ml
Gentamicine.....	0.10
Amphotericin.....	0.0005

Description /Technique

The purpose of transportation culture media is to maintain the viability of the virus contained in the sample, avoiding possible growth during transportation. Therefore, its formula is chemically defined and non-nutritive, to produce a reducing, buffered and osmotically balanced in the medium.

Proceed according to normative or methodology of the laboratory.

Quality control

Physical/Chemical control

Color : Light and transparent amber pH: 7.4 ± 0.2 at 25°C

Microbiological control

Not Performed - Only Sterility Test

Not Applicable

Microorganism

-

Growth

Not applicable

Sterility Control

Incubation 48 hours at 30-35 °C and 48 hours at 20-25 °C: NO GROWTH.
Check at 7 days after incubation in same conditions.

Bibliography

- HANKS, J.H. (1976) Hanks' Balanced Salt Solution and pH Control. Tissue Culture Association Manual. 3, 3. 2.
- HANKS, J.H. & R.E. WALLACE, (1949). Relation of Oxygen and Temperature in the Preservation of Tissues by Refrigeration. Proc.Soc. Exp. Biol. Med. 71, 196-200.
- RECHCIGL, M. Jr., (1977) CRC Handbook Series in Nutrition and Food. Section G: Diets, Culture Media, Foods Supplements. Volume IV. Culture Media for Cells, Organs and Embryos.
- CDC (2020) New Standard Operating Procedure for Creating Viral Transport Media