

## ● ● ● Autoclave Bags

**Technical Data Sheet:** Excelsior Order Code:  
ACB001E, ACB003E, ACB005E

### For Steam Autoclaves

#### Features

##### **Autoclave Bags**

- High Temperature Resistance
- Bags carry hazard warnings
- Autoclavable to 145°C
- Manufactured from heavy gauge polypropylene for superior performance.



#### Range of Sizes Available

<i>PRODUCT CODE</i>	<i>SIZE</i>	<i>BAGS PER BOX</i>	<i>DESCRIPTION</i>
ACB001E	310mm x 662mm	200	Small
ACB003E	413mm x 596mm	200	Medium
ACB005E	608mm x 816mm	200	Large

For additional product information:  
Please visit us at  
[www.excelsiorscientific.com](http://www.excelsiorscientific.com)  
Email us at [sales@excelsiorscientific.com](mailto:sales@excelsiorscientific.com)