



Puritan®
Quality since 1919

PurFlock Ultra® Flocked Swab

Ultrafine Tip with Polystyrene Handle - 100mm breakpoint



Description

An ideal choice in flock swabs with excellent absorption and elution of specimens when HydraFlock's high levels of absorbency aren't necessary. PurFlock Ultra's high purity characteristics makes it an ideal swab for various rapid diagnostic tests.

Patented: US, Canada, Europe, Australia, Japan, China
Patent pending: Hong Kong

Dimensions

6.02" Overall Length (152.908 mm)

Tip

- 0.140" Dia. (3.556 mm)
- 0.66" Length (16.76 mm)

Handle

- 0.1" Dia. (2.54 mm) tapered to 0.036" (.914 mm)
- 5.99" Length (152.146 mm)

Qualities

- Medical grade quality
- Highly absorbent tip
- Soft smooth tip
- Made in USA

Applications

- For nasopharyngeal and urethral cell collection
- Specimen collection
- Rapid diagnostic testing
- Forensic evidence collection



Ref #	Packaging
3317-U Non-Sterile	10/100 (10 bags per case, 100 per bag)
25-3317-U Sterile	10/50/1 (10 boxes of 50 individually wrapped swabs)