

# Sabouraud Dextrose Broth acc. to harm. EP/USP/JP

TN 1047 500 g



## Intended use

Test for *Candida albicans*. Culture medium for testing non-sterile products. This culture medium is based on the recommendations of the harmonised method acc. to harm. EP/USP/JP (2006).

## Form

Powder

## Farbe

Bright beige

## Storage

Dry, tightly closed, at 10...25 °C.

## Shelf life

3 years

## Directions

Suspend 30 g in 1 litre distilled water and autoclave for 15 minutes at 121 °C.

## Typical composition

Component	g/l
Pancreatic digest of casein	5.0
Peptic digest of animal tissue	5.0
Dextrose	20.0

## Final pH at 25 °C

5.6 ± 0.2

## Microbiological quality control

Test strain	Inoculum (CFU)	Growth	Incubation
<i>Candida albicans</i> ATCC 10231	10-100	good	aerobically, max. 72 h, 30-35 °C