

# Sabouraud Dextrose Agar acc. to harm. EP/USP/JP

TN 1264 500 g



## Intended use

Test for *Candida albicans*. Culture medium for testing non-sterile products. This culture medium is based on the recommendations of the harmonised method acc. to harm. EP/USP/JP (2006).

## Form

Powder

## Colour

Bright beige

## Storage

Dry, tightly, at 10...25 °C.

## Shelf life

5 years

## Directions

Suspend 65 g in 1 litre distilled water and heat until completely dissolved. Autoclave for 15 minutes at 121 °C.

## Typical composition

Component	g/l
Pancreatic digest of casein	5.0
Peptic digest of animal tissue	5.0
Dextrose	40.0
Agar	15.0

## Final pH at 25 °C

5.6 ± 0.2

## Microbiological quality control

Test strain	Growth	Colony morphology	Incubation
<i>Candida albicans</i> ATCC 10231	good	white, dry colonies	aerobically, max. 24 h, 30-35 °C