



## Product Specification Sheet

### **TS/15-A: NaCl**

Sterile 50cm<sup>2</sup> blue sponge dosed with 10 mlsaline solution – 0.9% Sodium Chloride (NaCl)in 100 ml plastic container

Product Code	TS/15-A: NaCl
Description	Sterile 50cm <sup>2</sup> blue sponge dosed with 10 mlsaline solution – 0.9% Sodium Chloride (NaCl)in 100 ml plastic container
Intended Use	Sampling of meat carcasses for microbiological analysis
Components	Cellulose blue sponge 25cm <sup>2</sup> x 2 sides ± 5mm. Blue plastic container with lid
Contents	0.9%NaCl - 10 ml Sodium Chloride Solution: De-ionised Water – 1 Litre NaCl–9g
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 50-92kGy
Shelf Life	12 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protectfrom light, heat and moisture in tightly closed original container
Pack Size	74 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.ii) Organism recovery (S.typhimurium)

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 3 Date : 27.02.2013

ProductSpec\_TS15-A+NaCl

