



## Product Specification Sheet

### **TS/15-T**

Sterile Blue Plastic Template with press 'hold' tab 10x10 cm

Product Code	TS/15-T
Description	Sterile Blue Plastic Template with press 'hold' tab 10x10 cm
Intended Use	Defined area Sample Collection. General laboratory use only
Components	Blue Plastic (soft) Medical Grade bag
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-46kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Precautions and Limitations	Template can only be used once or cross over contamination will occur
Pack Size	50 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 5 Date : 06.06.2016

ProductSpec\_TS15-T