



Product Specification - **TS/15-5B**.

Product Code	TS/15-5B
Description	Poultry boot swab (5 pairs) for litter analysis. Pre-Moistened in MRD (Maximum Recovery Diluent). Includes gloves & overboots in resealable outer bags.
Intended Use	Environmental sampling of poultry flocks for Salmonellae
Components	10x boot swabs – 2x Plastic overboots - 2x Plastic gloves - 2x resealable outer bags.
Contents	De-ionised water; Maximum Recovery Diluent*
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53KGy
Shelf Life	12 months from the date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Precautions & Limitations	Maximum Recovery Diluent has been modified to enhance enrichment of target organisms. Further selective steps will be required.
Pack Size	25 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility. ii) Organism Recovery

*adjusted to meet performance criteria

This specification is intended for guidance only. You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 0 Date: 16.04.2018
ProductSpec_TS15-5B