



TS/15-B:PSD

Blue Cellulose sponge in Easy Open Stomacher Pouch, individually packed & dosed with 10 ml of Peptone Saline Diluent (MRD) solution.

Product Code	TS/15-B:PSD
Description	Blue cellulose sponge in Easy Open Stomacher Pouch, individually packed & dosed with 10 ml of Peptone Saline Diluent (MRD) solution
Intended Use	Sampling of environmental areas and carcasses under special conditions (see Precautions & Limitations) for Microbiological analysis
Components	Cellulose blue sponge, Easy Open Stomacher Pouch
Contents	De-ionised Water – 1 Litre Maximum Recovery Diluent powder 9.5g
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	12 months from date of manufacture
Storage Conditions	4-25 °C; Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Precautions and Limitations	1. Some Fastidious organisms may struggle/fail to recover without extra supplement from the carcass/external sources.
Pack Size	200 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility. ii) Organism Recovery.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 4 Date : 06.10.2014
ProductSpec_TS15-B+PSD-V4

