



## Product Specification Sheet

### **TS/16-A**

Sterile dry swab with cotton mini tip and papershaft single in tube.

Product Code	TS/16-A
Description	Sterile dry swab with cotton mini tip and papershaft single in tube.
Intended Use	Culture swab; recommended for urethral sampling
Components	Swab with cotton mini tip (tip diameter 1.5 mm) and plasticised paper shaft, mounted in a polyethylene cap Polypropylene Tube Non-toxic glue Label & film
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-46kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Precautions and Limitations	
Pack Size	500 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 4 Date : 06.06.2016

ProductSpec\_TS16-A