



Product Specification Sheet

TS/3-A

Inoculating Needle (flexible), 20 in a sealed minigrip bag

Product Code	TS/3-A
Description	Inoculating Needle (flexible), 20 in a sealed minigrip bag
Intended Use	Sample throughout or Protect. General laboratory use only
Components	Plastic needle (flexible) Minigrip bag
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 20-46kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Precautions and Limitations	Laboratory use only
Pack Size	1000 in a box
Other	Quality control testing details on final product: Validation of irradiation sterility

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 5 Date : 06.06.2016

ProductSpec_TS3-A