



Product Specification Sheet

TS/30-B

L shaped Spreaders (firm), single in peel pouch

Product Code	TS/30-B
Description	L shaped Spreaders (firm), single in peel pouch
Intended Use	Sample throughout or Protect. General laboratory use only
Components	L shaped Plastic (firm) peel pouch
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-46kGy
Shelf Life	36 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Precautions and Limitations	Sterility is only guaranteed whilst bag is sealed. Once opened, any unused items should be discarded
Pack Size	1000 in a box
Other	Quality control testing details on final product: Validation of irradiation sterility

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 5 Date : 06.06.2016

ProductSpec_TS30-B