



## **Product Specification Sheet**

### **TS/4-E**

1 µl Yellow Culture Loops (soft), 50 in the sealed minigrip bag

Product Code	TS/4-E
Description	1 µl Yellow Culture Loops (soft), 50 in the sealed minigrip bag
Intended Use	Sample throughout. General laboratory use only
Components	1 µl Plastic Yellow Loops (soft) Minigrip bag
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 20-46kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Precautions and Limitations	Laboratory use only
Pack Size	1000 in a box
Other	Quality control testing details on final product: Validation of irradiation sterility

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 5 Date : 06.06.2016

ProductSpec\_TS4-E