



## Product Specification Sheet

### **TS/5-10**

5ml Amies with charcoal transport medium in tube (blue cap) with orange cap cotton mini tipped plasticised paper shaft swab in a peel pouch

Product Code	TS/5-10
Description	5ml Amies with charcoal transport medium in tube (blue cap) with orange cap cotton mini tipped plasticised paper shaft swab in a peel pouch
Intended Use	Bacteriological samples only. Collection and transport of clinical specimens
Components	Blue cap; Polypropylene Tube; Orange cap; cotton tip plasticised paper shaft swab; Non-toxic glue; Label & film; Amies Transport Medium;
Contents	De-ionised water; Amies Transport Medium;
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-46kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Precautions and Limitations	Recovery after 48 hours will be unpredictable and target organisms will be lost. Fatty acid sensitive strains will give decreased recovery or be stressed in the presence of cotton
Pack Size	100 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.ii) Organism Recovery

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 4 Date : 06.06.2016

ProductSpec\_TS5-10