



Product Specification Sheet

TS/6-S6

Sterile 6" sharpened tip cotton swab with a wooden shaft in 6" polypropylene outer tube

Product Code	TS/6-S6
Description	Sterile 6" sharpened tip cotton swab with a wooden shaft in 6" polypropylene outer tube
Intended Use	Culture swab; general laboratory use; forensic or precision work
Components	6" Polypropylene Tube 6" Swab with Blue cap Cotton tip/wooden shaft Non-toxic glue Label & film
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-46kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Precautions and Limitations	
Pack Size	500 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 3 Date : 06.06.2016

ProductSpec_TS6-S6