



Product Specification Sheet

TS/9-B

Sterile swab with a wooden shaft and cotton tip, 5 in the peel pouch.

| | |
|-----------------------------|---|
| Product Code | TS/9-B |
| Description | Sterile swab with a wooden shaft and cotton tip, 5 in the peel pouch. |
| Intended Use | Culture swab & general laboratory use |
| Components | Cotton tip/wooden shaft swab Peel pouch Non-toxic glue Label & film |
| Contents | Not applicable. |
| Method of Sterilisation | Sterilised by gamma irradiation at dosage of 25-46kGy |
| Shelf Life | 24 months from date of manufacture |
| Storage Conditions | 4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container |
| Precautions and Limitations | Once opened, pack must be used aseptically & discard any swabs left over |
| Pack Size | 1000 in a box |
| Other | Quality control testing details on final product:i) Validation of irradiation sterility. |

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 5 Date : 06.06.2016

ProductSpec_TS9-B