



Product Specification Sheet

TS/70-H

Aluminium Cryoblock for optimal PROTECT use

Product Code	TS/70-H
Description	Aluminium Cryoblock for optimal PROTECT use
Intended Use	Maintaining longevity of microorganisms held in protect vials whilst in use (ie outside the freezer).
Components	Aluminium block, polyurethane insulation, Non-toxic glue
Contents	Aluminium block, polyurethane insulation, Medium Minigrip bag 203mm x 305mm
Method of Sterilisation	N/A
Shelf Life	N/A
Storage Conditions	N/A
Precautions and Limitations	
Pack Size	1 in minigrip bag
Other	Quality control testing details on final product.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 2 Date : 06.06.2016

ProductSpec_TS/70-H