

HAND DISINFECTION



1-LITRE TRIGGER SPRAY HANDS-FREE DISPENSER

DESCRIPTION

Stainless steel AISI 304 BA hands-free dispenser. Designed for wall-mounting and fitting most of the European 1-litre trigger sprays on the market. The elongated lever allows operators to disinfect their gloved hands, by using the arm or the elbow, avoiding any cross-contamination. Easily sanitizable and suitable for autoclaving.

HAND DISINFECTION

BENEFITS

CROSS-CONTAMINATION PREVENTION

The design of this support allows the dispensing of the disinfectant without the hands' aid, to prevent the risk of cross-contamination.

CONVENIENT

The system allows to dose the right quantity of product for gloved hands disinfection, without any waste.

EASY TO USE

The dispenser, in polished stainless steel, is installed on the wall and can be easily removed for autoclaving. The design allows easy cleaning.

PRODUCT		
CODE	DESCRIPTION	PACKAGING
AM820	1-LITRE TRIGGER SPRAY HANDS-FREE DISPENSER	1 piece

COMPANY WITH
QUALITY SYSTEM
CERTIFIED BY DNV GL
= ISO 9001:2015 =

AM Instruments srl

Via Isonzo, 1/c 20812 Limbiate (MB) - Italy
Tel. +39 02 872892.1 - Fax +39 02 872892.00
aminstruments.com

