

Product specification MSC93050 – Oral Flocked swab

Product description

MSC-93050 is a flocked swab of sampling for Laryngopharyngeal clinical diagnostic, it utilizes state of the art "spray on technology" that the flocking process by means of an electro-static charge perpendicularly attaches millions of nylon microfibers on the medical grade handle tip. The flocked swab is ideal for collecting large amount of cells and rapid elution of the specimens that instantly releases the cells into the transport medium. It has been well recognized and adopted by the diagnostic test kit manufacturers who produce reagents in molecular genetics, forensics and clinical laboratories sectors.



Product Specifications

Item Number	Total length	Breakpoint	Tip Thickness	Tip Width	Tip length	Handle Diameter
MSC-93050	155mm	81mm	5mm	5mm	16mm	2.5mm

MSC-93050 Flocked swabs utilize innovative nylon fibre technology to provide outstanding sample collection and release. Nylon fiber swab that allows a better clinical sample collection and the following immediate release into the liquid transport medium. It has a plastic shaft with "breakpoint" to allow the swab break into the test tube and successive collection by screw "capture" cap after the tube closing.

Sterilization Process

We take ETO sterilization for MSC-93050.