



**TS/15-B:NaCl**

Blue Cellulose sponge in Easy Open Stomacher Pouch, individually packed & dosed with 10 ml of 0.9% saline solution.

Product Code	TS/15-B:NaCl
Description	Blue cellulose sponge in Easy Open Stomacher Pouch, individually packed & dosed with 10 ml of 0.9% saline solution
Intended Use	Sampling of Meat Carcasses for Microbiological analysis
Components	Cellulose blue sponge, Easy Open Stomacher Pouch
Contents	De-ionised Water – 1 Litre Sodium Chloride 9.0g
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	12 months from date of manufacture
Storage Conditions	4-25 °C; Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Precautions and Limitations	1. Some Fastidious organisms may struggle/fail to recover without extra supplement from the carcass/external sources.
Pack Size	200 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility. ii) Organism Recovery.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 3 Date : 06.10.2014

ProductSpec\_TS15-B+NaCl\_V3

