



Product Specification Sheet

TS/15-CP

Foam tipped swab with a blue flexible polypropylene shaft, in polypropylene outertube.

Product Code	TS/15-CP
Description	Foam tipped swab with a blue flexible polypropylene shaft, in polypropylene outertube.
Intended Use	Environmental and drain sampling.
Components	Jumbo blue cap 12" Straw Foam Sponge Head 12" Plugged Tube Dimensions: Shaft 295+5 mm x 16-18 mm at widest point (swab head diameter)
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-46kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Precautions and Limitations	Do not use in solvents
Pack Size	100 in a box
Other	Quality control testing details on final product: Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 4 Date : 06.06.2016

ProductSpec_TS15-CP