



Product Specification Sheet

TS/16-S

Sterile dry swab, single in peel pouch. Cotton mini tip and paper shaft swab, breakpoint at 45mm.

Product Code	TS/16-S
Description	Sterile dry swab, single in peel pouch. Cotton mini tip and paper shaft swab, breakpoint at 45mm.
Intended Use	Culture swab; recommended for urethral sampling
Components	Swab with cotton mini tip (tip diameter 1.5 mm) and plasticised paper shaft, breakpoint at 45mm Plastic peel pouch Non-toxic glue Label &film
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-46kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Precautions and Limitations	Recovery after 24 hours will be unpredictable and target organisms will be lost.
Pack Size	1000 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 5 Date : 06.06.2016
ProductSpec_TS16-S