



## Product Specification Sheet

### **TS/19-G**

Swab with a breakable polystyrene shaft and polyester tip in a polypropylene tube

Product Code	TS/19-G
Description	Swab with a breakable polystyrene shaft and polyester tip in a polypropylene tube
Intended Use	Culture swab; Patient sampling & general Laboratory use
Components	Polypropylene Tube Cap polyester tip/polystyrene shaft swab Non-toxic glue Tube label Box label & film
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-46kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Precautions and Limitations	
Pack Size	500 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 3 Date : 06.06.2016

ProductSpec\_TS19-G