



## Product Specification Sheet

### **TS/4-C**

1 µl Clear Culture Loops (firm), single in pouch peel

Product Code	TS/4-C
Description	1 µl Clear Culture Loops (firm), single in pouch peel
Intended Use	Sample throughout. General laboratory use only
Components	1 µl Plastic Clear Loops (firm) peel pouch
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-46kGy
Shelf Life	36 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Precautions and Limitations	Laboratory use only
Pack Size	1000 in a box
Other	Quality control testing details on final product: Validation of irradiation sterility

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 5 Date : 06.06.2016

ProductSpec\_TS4-C