



Product Specification Sheet

TS/5-18

5ml Amies with charcoal transport medium in tube with black cap pure viscose tipped (polystyreneshaft) swab in a peel pouch Latex free

Product Code	TS/5-18
Description	5ml Amies with charcoal transport medium in tube with black cap pure viscose tipped (polystyreneshaft) swab in a peel pouch Latex free
Intended Use	Bacteriological samples only. Collection and transport of clinical specimens
Components	Black cap Polypropylene Tube Black cap viscose/polystyrene swab Non-toxic glue Label & film Transport medium
Contents	De-ionised waterAmies Charcoal Transport MediumpH – 7.15 ± 0.2
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-46kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Precautions and Limitations	Recovery after 48 hours will be unpredictable and target organisms will be lost. Fatty acid sensitive strains will give decreased recovery or be stressed in the presence of cotton
Pack Size	100 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.ii) Organism Recovery

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.