



Product Specification Sheet

TS/5-18A

4ml Amies with charcoal transport medium in tube with black cap pure cotton tipped polystyrene shaftswab in a peel pouch

Product Code	TS/5-18A
Description	4ml Amies with charcoal transport medium in tube with black cap pure cotton tipped polystyrene shaftswab in a peel pouch
Intended Use	Bacteriological samples only. Collection and transport of clinical specimens
Components	Black cap Polypropylene Tube (80 mm) Black cap cotton tipped plastic swab (shaft length 80 mm) Non-toxic glue Label & film
Contents	De-ionised water Amies Charcoal Transport Medium pH – 7.15 ± 0.2
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-46kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Precautions and Limitations	Recovery after 48 hours will be unpredictable and target organisms will be lost. Fatty acid sensitive strains will give decreased recovery or be stressed in the presence of cotton
Pack Size	100 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility. ii) Organism Recovery (N.gonorrhoeae, H.influenzae)

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.