



## Product Specification Sheet

### **TS/6-A**

Sterile swab with a wooden shaft and cotton tip, in polypropylene outer tube

Product Code	TS/6-A
Description	Sterile swab with a wooden shaft and cotton tip, in polypropylene outer tube
Intended Use	Culture swab; general laboratory use
Components	Polypropylene TubeSwab with Blue cap Cotton tip/wooden shaft Non-toxic glue Label & film
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-46kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Precautions and Limitations	Unsuitable for certain sample areas due to wood shaft and possibility of splinters.
Pack Size	500 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 5 Date : 06.06.2016

ProductSpec\_TS6-A