



Product Specification Sheet

TS/133-M

Autoclave Bag for high temperature (<150 °C) autoclaving printed with biohazard symbols

Product Code	TS/133-M
Description	Autoclave Bag for high temperature (<150 °C) autoclaving printed with biohazard symbols
Intended Use	High temperature (<150°C) autoclaving of hazardous material
Components	Flexible clear polypropylene (PP) film with welded sides, 38 micron thick, 800mm x 600mm (± 10mm) with blue Hazard print on one side
Contents	Not applicable as supplied empty.
Method of Sterilisation	N/A - supplied unsterile
Shelf Life	N/A
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	500 in a box
Other	

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 10.02.2014

ProductSpec_TS133-M



Intertek



Intertek



UKAS
QUALITY
MANAGEMENT

061